

Quiz: Reprocessing of Medical Equipment

Name _____

1. What information is important to read when looking at a disinfectant label?
 - A. Manufacturer's instructions for use
 - B. Contact time to follow
 - C. Type of germs the disinfectant can kill, e.g., bacteria, virus, spores
 - D. All of the above
2. Transfer lifts are considered noncritical equipment and should be disinfected with?
 - A. Low level disinfectant
 - B. Intermediate level disinfectant
 - C. High level disinfectant
 - D. Only needs cleaning not disinfecting
3. Which of the following is true regarding cleaning and disinfecting reusable medical equipment?
 - A. Cleaning must always be done before disinfecting
 - B. Reprocessed medical equipment and soiled equipment can be stored together
 - C. Keep the surface of medical equipment wet with a disinfectant for 30 seconds
 - D. None of the above
4. True or False: Nursing staff only need to wear personal protective equipment (PPE) if cleaning and disinfecting reusable medical equipment of a resident on transmission-based precautions
5. True or False: Environmental Services staff may be assigned to clean and disinfect reusable medical equipment.